

Autoclavable Syringe Components

Button Valves

23.1232.01 Buttons (2)



- Requires spring 013.084.01 (10)

2-Button Conversion Kit

23.1194.00



- Converts from single autoclavable button to two autoclavable buttons
- Includes: (2) valve assemblies, (1) spacer, screw, (2) springs, silicone lubricant, and hex key

Service Tips

USE ONLY A-DEC SILICONE LUBRICANT (98.0090.01)

Using petroleum-based lubricants (jellies, creams, or lotions) causes the O-rings in the syringe to swell and/or deteriorate, resulting in syringe malfunction or failure.

USE ONLY A-DEC SYRINGE TIPS

Using disposable syringe tips or tips manufactured by a company other than A-dec may result in damage to the O-rings inside the tip retainer nut assembly. Damaged O-rings may result in syringe tip ejection.

Identify genuine A-dec tips by the "a" located above three grooves on each tip.

Autoclavable Syringe Valve Kit

23.1234.00



- For syringes beginning October 1995
- Includes button valves, springs, retaining screw, hex key, and silicone lubricant

Autoclavable QD Syringe Service Kit

23.1235.01



- Includes syringe tip, screw, pin, syringe nut with spacer, springs, button valves, stem assemblies, terminal barbs, silicone lubricant and O-rings

Button Valves

23.1021.01 Buttons (2)



- Requires spring 013.084.01 (10)