Handpiece Holder Valves

2-Way Air Bleed Valve for Cascade Unitized Holder

33.0132.00 2-way valve 99.0590.00 Actuator only



- · Use for air bleed
- . Snap fit mounting to actuator
- 1/8" barb
- · Actuator only



. SERVICE TRANS

2-Way Air Bleed Valve for Excellence Holder

33.0025.00



- · Use for air bleed
- · Snap fit mounting
- 1/8" barb

2-Way Air Bleed Valve

33.0025.01



- · For Continental whip holder
- · Use for air bleed
- · Snap fit mounting
- 1/8" barb

2-Way Air Bleed Valve for Cascade Individual Holders

33.0025.00 2-way valve

99.0590.00 Actuator only

99.0592.00 Actuator w/valve





- · Use for air bleed
- . Snap fit mounting to actuator
- · 1/8" barb

儒

Traditional and Continental

99.0654.02 Actuator, Surf 4


99.0654.03 Actuator, White 2

Handpiece Actuator Valve and Tubings

 For A-dec delivery system models: 532, 533, 542, 332, 333, 334, and 335 handpiece holders and whites

2- or 3-Way Micro Valves

33.0134.00 2-Way (brass) 33.0138.00 3-Way (SST)



- · 2-Way, brass ball
- 3-Way, stainless steel ball (SST), replaces 29.0185.00
- 1/8" barb

Actuator Valve, 12 O'Clock Doctor's Holder & Tubing

99.0696.01 541 Actuator, Surf 4 99.1135.01 541 Actuator, White 2



 For Doctor's handpiece holders on a 541 12 o'clock delivery system

Performer Holder Valve & Tubing

99.0724.00



- · 3-Way micro valve with tubing
- Retrofits to replace previous holder valves
- Replace valve by detaching tubing from control block

Note

Micro valve 29.0185.00 has been replaced by 3-way micro valve 33.0138.00 (see middle column, middle product).