Snap Ring Pliers

009.007.00



- Simplifies installation and removal of both internal and external snap rings
- Fits all snap ring sizes used in A-dec equipment

Umbilical Tubing Stringer

009.015.00



- Use to string new or additional tubing and/or wiring into an existing umbilical assembly
- 12' of stringer with threading holes at each end

Tubing Pliers

009.014.00



 Simplifies installation of 1/8" OD tubing to barbed fittings

A-dec Sleeve Tool

98.0072.00



 Use to secure the 1/4" and 1/8" tubing sleeves when installing tubing on barb fittings

O-Ring Tool Kit

009.013.00



- · For quick field repair
- Includes three tools, which fit the four smallest 0-ring sizes used in A-dec equipment

Tubing Cutters

009.036.00



Hemostat

009.008.00



 Temporarily stops air or water flow through tubing while troubleshooting or repairing the unit

Syringe Wrench

98.0900.00



Tubing Strippers

009.035.00

