
# Century Plus

# Century Plus Control Block

38.0540.00



· With pigtails

### Century Plus Cartridges

38.0517.00 Air bleed 38.0518.00 Check valve 38.0520.00 Water valve



Air bleed (brass base)



Check valve (blue base with duckbill)



(red base)

# Drive Air Flow Adjustment Stem

38.0510.00 Stem 035.034.01 0-ring (10)



· Includes O-ring

# Water Coolant Flow Adjustment Stem


38.0516.00 Stem 035.034.01 0-ring (10)



· Includes 0-ring

#### Syringe Flow Control Barb

38.0555.00



- · For air and water
- · 1/8" barb
- · 10-32 threaded fitting

#### Air Coolant Flow Adjustment Stem

38.0526.00 Stem 030.003.02 0-ring (10)

· Includes 0-ring



#### Handpiece Control Adjustment Keys

98.0867.02 (5)



- Autoclavable adjustment keys
- For use in Cascade and Decade 4400 handpiece controls

#### Flush Valve Stem

38.0525.00 Stem

035.033.01 O-ring (10)

034.001.01 O-ring (10)

013.072.00 Spring



# Gauge, Horizontal

026.083.01



- Panel mount
- 0 60 psi
- · 1/2" x 13/4"
- · 6-32 mounting holes
- 10-32 port
- · Includes 1/8" barb